J. Alourong Aren um On

PEIROTECH PEIROTECH PEIROTECH

## FIELD REPORT

TO

:

J. Hagan.

**FROM** 

R. Rosas / G. Ruesta

**SUBJECT** 

Well: LO16-14 Lobitos Area

Treatment with Bacteria (PARA BAC &

CORROSO BAC).

DATE	TIME	DESCRIPTION	
Abr.,20, 99	14:00	Arrived to LO16 Platform.	
: : :	14:00 - 19:00	Hook up Halliburton lines and tested with 3,800 psi, on well LO16-14.	
1. tim to		Opened well: 360 psi. Pumped as follows: - 80 bbls. of fresh water w/ 500 lbs. KCL mixed with:	
rulation to op. Banal Sa Sn isolated	zline	- 03 drums of PARA BAC-S 01 drum of PARA-BAC-X.	
TV inlated		- 01 drum of CORROSO BAC. Displaced with 63 Bbls. of Diesel.	
Sn 1800	/	Pump Pressure: 60 psi. Average Rate: 3 BPM. Final Pressure: 1,200 psi.	14 1 2 <i>6</i>
		Final Pressure: 1,200 psi.  Shut in well ( will continue for 5 days ).— it will open on Rigged down Halliburton lines and towed to Talara.	at 19 hrs.
i		Note: Pressure began to rise after pumped 68 bbls.	

cc: J. Meyers

R. Francia

O. Nieto

Well File

5